

CODAN CYTO[®] Product Line



CODAN CYTO[®]

A Closed Solution for Cytotoxic Drug Delivery that is Clean, Safe and Simple

CODAN's Closed Solution calls for:

1. Preparation with CODAN Vial Spikes under a Class III enclosed, ventilated hood (BSC)
2. Closed-system administration using neutral solution barriers and no-disconnect products


Vial Spikes	Primary Connect Sets	Secondary Connect Sets	IV Extension Sets	Gravity Sets	Accessories
C-31C	CYTO-16A	CYTO-21V	CYTO-1BC	CYTO-13G	BC-128
C-32C	CYTO-17A	CYTO-22V	CYTO-2BC	CYTO-14G	BC-129
C-34C	CYTO-18A	CYTO-23V	CYTO-3BC	CYTO-41T	BC-130
C-350			CYTO-5BC		
C-355					

Vial Spikes:

Featuring 0.2 micron air eliminating filter with 99.999% aerosol retention

C-350


25.7280



Universal Vial Spike with 0.2 micron air eliminating filter, female luer-lock with cap, 50/cs

C-355

16.7102



Chemoprotect[®] Spike with 0.2 micron air eliminating filter, 5.0 micron fluid filter, female luer-lock with cap, 50/cs

C-31C

25.7141



Universal Vial Spike with swabable needle-free adapter, 0.2 micron air eliminating filter, 50/cs

C-32C

25.7321



Chemoprotect[®] MicroSpike with 0.2 micron air eliminating filter, 5.0 micron fluid filter, female luer-lock with cap, 50/cs

C-34C

25.7242




Chemoprotect[®] MicroSpike with 0.2 micron air eliminating filter, 5.0 micron fluid filter with swabable needle-free adapter, 50/cs

CODAN CYTO[®] - From Pharmacy to Chair Side

Primary Connect Sets:


CYTO-16A

76.7221 Primary Connect Set with universal spike adapter, for connection to any primary pump set in the market, with 2 check valves as access ports and slide clamp, 50/cs




CYTO-17A

76.7241 Primary Connect Set with universal spike adapter, for connection to any primary pump set in the market, with 4 check valves as access ports and slide clamp, 50/cs



CYTO-18A

76.7261 Adapter Set with non-vented spike, slide clamp, inline spike adapter and Y-site check valve as access port, 50/cs



Secondary Connect Sets:

CYTO-21V

76.7042 Secondary Connect Set with luer activated valve, slide clamp and CODAN OnceLock with FlowStop Cap, for connection to Primary Connect Set's access port, Chemoprotect[®] Transport & Waste Bag included, 50/cs



CYTO-22V

76.7041 Secondary Connect Set with luer activated valve, slide clamp, 0.2 micron air eliminating filter and CODAN OnceLock with FlowStop Cap, for connection to Primary Connect Set's access port, Chemoprotect[®] Transport & Waste Bag included, 50/cs



CYTO-23V


76.7182 Secondary Connect Set with luer activated valve, inline chamber, slide clamps and CODAN OnceLock with FlowStop Cap, for connection to Primary Connect Set's access port, Chemoprotect[®] Transport & Waste Bag included, 50/cs



IV Extension Sets:

CYTO-1BC


76.7044



Safety IV Extension Set with StopPost-StopCock, SpinLock with FlowStop Cap and CombiLock DeadEndEr, 50/cs

CYTO-2BC


72.7065



Bifurcated IV Extension Set with 2 check valve needle-free access ports, 0.2 micron air eliminating filter, needle-free Y-connector, pinch clamps and SpinLock with FlowStop Cap, 50/cs

CYTO-3BC


76.7082



Bifurcated IV Extension Set with male luer-lock and 2 female check valves as access ports, 50/cs

CYTO-5BC

72.7681




IV Extension Set with check valve needle-free access port, slide clamp, 0.2 micron air eliminating filter, needle-free Y-injection site and spin-lock with FlowStop Cap, 50/cs

Gravity Sets:

CYTO-13G


76.7083



Gravity Set with inline check valves, 20 drop inline chamber with 15 micron filter, slide clamp, roller clamp, pinch clamp, 0.2 micron air eliminating filter, needle-free Y-injection site, spin-lock with FlowStop Cap and 2 check valves as access ports, 50/cs

CYTO-14G


76.7093



Gravity Set with inline check valves, 20 drop inline chamber with 15 micron filter, slide clamp, roller clamp, pinch clamp, 0.2 micron air eliminating filter, needle-free Y-injection site, spin-lock with FlowStop Cap and 4 check valves as access ports, 50/cs

CYTO-41T

76.7034




Gravity Set with 20 drop inline chamber with 15 micron filter, pinch clamps, roller clamp, 0.2 micron air eliminating filter, needle-free Y-injection site, spin-lock with FlowStop Cap and one Y check valve as access port, 50/cs

Accessories:

BC-128


16.7253



FlowStop Cap, allows easy priming when attached to various male luer-locks at the bottom of a pump set, 504/cs

BC-129

16.7254



CYTO Clip 2, for connection between infusion bags and spike of secondary connect sets, to prevent tampering and as a reminder of the CODAN "No Disconnect" philosophy. Compatible with Baxter Viaflo (new version), Fresenius Kabi, Maco Pharma N and PAB® Container B. Braun USA. 100/dispenser box, 500/cs

BC-130

16.7255




CYTO Clip 3, for connection between infusion bags and spike of secondary connect sets, to prevent tampering and as a reminder of the CODAN "No Disconnect" philosophy. Compatible with Baxter Viaflex® and Hospira USA. 100/dispenser box, 500/cs

PAB® is a registered trademark of B. Braun Medical, Inc.
Viaflex® is a registered trademark of Baxter International, Inc.



CODAN CYTO® HighLights:

- A Closed Solution that is Clean, Safe and Simple
- Engineered Controls
- Universal Spike Adapter
Compatible with all Pump Systems
- System Greatly Reduces Sympathetic Flow
- No Disconnects = No Exposure
- No Back Priming
- No Accidental Needle Sticks
- FDA, CGMP, ISO, CE and CMDCAS Standards



For additional information about **CODAN CYTO®** and other CODAN products, call us toll-free at **(800) 33-CODAN (26326)** or visit our website at **www.codanusc corp.com**.